



2024.10.30

Vision-Based Intelligent Gait Environment Detection Algorithm for Assisting Stroke Patients with a Soft Ankle Exosuit

Gayoung Yang and Brian Byunghyun Kang

Intelligence Robotics Lab

School of Intelligent Mechatronics Engineering

Sejong University



세종대학교
SEJONG UNIVERSITY



Intelligent Robotics
Laboratory

I. Introduction

II. Soft Exosuit

1. Soft Exosuit Prototype
2. Control Strategy

III. Environment Recognition System

IV. Experiments

V. Results


1. Real-Time Environment Classification Performance
2. Joint Kinematics
3. Muscle Reduction

VI. Conclusion

I. Introduction


Stroke Patients

- functional impairments & gait disturbances
- **foot drop** : reduction in foot clearance (swing phase of gait cycle)



stair ambulation

- common activity
- challenging task : joint range, muscular moment $\uparrow$
- ascent: insufficient clearance between foot & stair nosing



Lower Limb Exoskeleton

- Assisting gait and rehabilitation
- Crucial for activities of daily living (ADL)

I. Introduction

Exoskeleton

inertial measurement unit (IMU) & insole sensor

- provide real-time gait phase information
- Limitation: Lack of prior environmental information

vision – based perception



human
peripheral vision



exoskeleton
vision system

- ✓ visual recognition of step location
- ✓ continuous monitoring stage

- ✓ awareness of surrounding environment
- ✓ locomotion guidance

I. Introduction

Deep learning-based methodologies

influence robotics field

environment classification + exoskeleton

[challenge]

- transition phase accuracy
- stroke patient consideration




Vision-based Intelligent gait environment detection algorithm for ankle exosuit

- Integration of deep learning-based environment classification into exoskeleton application
- detailed gait environment classification (YOLOv8) + spatio-temporal feature extraction (LSTM)
- soft ankle exosuit aimed at stroke patients


II . Soft Exosuit

Exosuit Design



- flexible and lightweight materials
- Bowden cable-driven method

Control Strategy



flowchart of gait-environment-based exoskeleton control strategy

III . Environment Recognition System

Data Acquisition · Class Definition

[customized dataset]

- level ground and stair environment
- 960 x 540 pixel
- training: 16,000 | valid: 2,000 | test: 2000

-> 8:1:1

train

valid

test



IC prev


initial incline

IC

mid-incline

IC next

final incline



LG

level ground



DC prev

initial decline

DC

mid-decline

DC next


final decline

III . Environment Recognition System

YOLOv8 Environment Classification

- object detection framework: one-stage detection approach
- provide real-time contextual insight for frame

{ Rapid Detection: critical for detecting gait environment transitions
Spatial analysis: Effective in analyzing global spatial information



III . Environment Recognition System

LSTM-Based Spatio-Temporal Feature Extraction

- extract temporal features: adapting to environment transitions



Data Preparation

- Frame Sampling: **20 FPS**
- 1 class for 10 frames,
with a stride of 5 (**4 labels/sec**)
- Sequence shape: **(1714, 10,8)**




Class Labeling

- **Most probable Class:**
Chosen from YOLO frame-wise probabilities
- **Maximum Voting & Average Voting**




Model Training

- 100 epochs
- batch sizes : 32
- learning rate: 0.001
- Adam optimizer



IV . Experiments


Experimental Setup



Condition





- 1) without exosuit (*No Exo*)
- 2) Exosuit worn but unpowered (*Exo Off*)
- 3) Exosuit worn and powered (*Exo On*)

10-minute familiarization with exosuit
20-minute rest between conditions



V . Results

Real-Time Environment Classification Performance



YOLO

Low Recall in Transition Classes (IC_{prev} , DC_{prev})

Overall F1 score : **0.774** - transition misclassification

YOLO - LSTM

Overall F1 score : **0.971** - balanced classification performance

YOLO




Misclassified LG as stairs 1.5s (6 frames) before stair transition

YOLO - LSTM

rapid and accurate classification during the transition phase

V . Results

Joint Kinematics



peak ankle joint dorsiflexion angle

	LG	IC	DC
<u>no exo</u>	8.88° ± 3.95°	7.92° ± 1.48°	9.14° ± 2.27°
<u>exo off</u>	6.74° ± 2.08°	25.99° ± 2.78°	25.08° ± 3.63°
<u>exo on</u>	27.13° ± 5.71°	21.06° ± 2.08°	25.65° ± 5.16°




LG: **significant increase** with the *exo on*

IC: **increase** with the *exo on*, *exo off*

V . Results

Muscle Reduction






average TA muscle activation

	LG	IC	DC
<u>no exo</u>	20.63 ± 11.56%	18.44% ± 9.92%	10.75% ± 4.25%
<u>exo off</u>	12.77% ± 5.02%	15.24% ± 12.77%	12.34% ± 10.93%
<u>exo on</u>	13.22% ± 3.43%	9.79% ± 10.12%	12.09% ± 4.81%

LG, IC: **reduction** in *exo off*, *exo on* condition
IC: slight increases with exosuit

V . Results

Muscle Reduction



average MG muscle activation

	LG	IC	DC
<u>no exo</u>	7.88% ± 6.61%	11.71% ± 8.83%	9.97% ± 7.35%
<u>exo off</u>	8.51% ± 4.04%	10.66% ± 5.99%	7.62% ± 3.41%
<u>exo on</u>	4.87% ± 1.21%	4.70% ± 1.69%	6.41% ± 2.19%

LG, IC, DC: **reduction** in *exo on* condition

Conclusion

- **vision-based detection algorithm** for an ankle exosuit, focusing on improving gait support for stroke patients
- integration of **YOLOv8** for environment classification and **LSTM** for spatio-temporal recognition
- **rapid and accurate classification** during gait environment transition
- **increased ankle dorsiflexion / Reduced TA and MG muscle activation**

Future work

- Validate across diverse real-world scenarios and multiple stroke patients
- Enhance the effectiveness and generalization ability of the exosuit by improving environment classification algorithm



actuators

an Open Access Journal by MDPI



Special Issue: Recent Advances in Soft Actuators, Robotics and Intelligence



Contact us:

Prof. Dr. Brian Byunghyun Kang

brianbkang@sejong.ac.kr

Actuators Editorial Office actuators@mdpi.com

Deadline for manuscript submissions: 31 March 2025

https://www.mdpi.com/journal/actuators/special_issues/L21IU2ZFN9

Call for Paper

This Special Issue will focus on the latest advances in the field of **soft technologies for robot realization**. Here, we invite you to share your latest findings by submitting high-quality manuscripts that address the following challenges in the field:

- Soft actuators and pliable mechanisms in robotics;
- Artificial muscles;
- Material and structural design of soft actuators;
- Soft materials in soft actuators or robots;
- ...



2024.10.30

Thank You

Gayoung Yang

Intelligence Robotics Lab

School of Intelligent Mechatronics Engineering

Sejong University



세종대학교
SEJONG UNIVERSITY



Intelligent Robotics
Laboratory